

ROTANODE™  
E7100X  0197

**Rotating Anode X-ray Tube Assembly**

- ◆ High speed rotating anode X-ray tube assembly for high energy radiographic operations.
- ◆ For the purpose of general diagnostic X-ray procedures.
- ◆ This tube has foci 1.2 and 0.6, and is available for a maximum tube voltage 150 kV.
- ◆ This tube housing assembly has specially processed rhenium-tungsten faced molybdenum target of 100 mm diameter anode disc and is accommodated with IEC60526 high-voltage cable receptacles.



**General Data**

**IEC Classification (IEC60601-1:2005) ..... Class I ME EQUIPMENT**

**Electrical:**

Circuit:

High Voltage Generator ..... Constant Potential High-Voltage Generator  
Grounding ..... Center-Grounded

Nominal X-ray Tube Voltage (IEC60613:2010):

Radiographic ..... 150 kV  
Fluoroscopic ..... 125 kV

Nominal Focal Spot Value (IEC60336:2005):

Large Focus ..... 1.2  
Small Focus ..... 0.6

Nominal Anode Input Power (at 0.1s) ..... See rating charts

		50 Hz	60 Hz	180 Hz
Large Focus	.....	55.5 kW	59 kW	100 kW
Small Focus	.....	21.5 kW	24 kW	40 kW

Nominal Radiographic Anode Input Power (IEC60613:2010):

		50 Hz	60 Hz	180 Hz
Large Focus	.....	44 kW	48 kW	80 kW
Small Focus	.....	17 kW	18 kW	30 kW

★The information contained herein is presented only as a guide for the applications of our products.  
No Responsibility is assumed by TOSHIBA ELECTRON TUBES & DEVICES CO.,LTD.(TETD) for any infringements of patents or other rights of the third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of TETD or others.  
★The information contained herein may be changed without prior notice. It is therefore advisable to contact TETD before proceeding with the design of equipment incorporating this product.

Motor Ratings:<sup>1)</sup>

Stator: XS-AG

	Starting		Running	
	180	50/60	180	50/60
Driven Frequency [Hz]	180	50/60	180	50/60
Input Power [W]	3600	1500	200	100
Voltage <sup>4)6)</sup> [V]	370	170	95	50
Current <sup>5)</sup> [A]	10.3	10.6	2.7	2.6
Min. Speed Up <sup>2)8)</sup> [s]	1.2	0.8	-	-
Capacitor [ $\mu$ F]	6	44	6	44
Min. Braking <sup>3)8)</sup> [s]	3.0 (DC 80V)			

Note 1) To be obtained with TOSHIBA starter RS-200 or equivalent.

2) The speed up time from normal speed to high speed is 2/3 times of the specified speed up time from 0 to high speed, which is described on motor rating table.

3) To be applied for high speed rotation.

4) Applied voltage between common and main terminal.

5) Common current.

6) The every applied voltage must be never exceeded 110% of the above specification.

7) No more than two high speed starts per minute are permissible.

8) The speed-up time is allowed up to 110% of the above specification.

## Anode Speed:

50 Hz .....	Minimum 2700 min <sup>-1</sup>
60 Hz .....	Minimum 3200 min <sup>-1</sup>
180 Hz .....	Minimum 9700 min <sup>-1</sup>

## Stator Resistance:

Common-Main Winding .....	9.4 $\Omega$
Common-Auxiliary Winding .....	28.3 $\Omega$
Resistance between Housing and Low Voltage Terminals .....	Minimum 2 M $\Omega$
Normal Operating Range of the Housing Temperature .....	16 ~ 75 °C
Mode of Operation .....	Intermittent

**Mechanical:**

Dimensions .....	See dimensional outline
Overall Length .....	496 mm
Maximum Diameter .....	195 mm
Target:	
Anode Angle .....	12 degrees
Diameter .....	100 mm
Construction .....	Rhenium-Tungsten faced Molybdenum
Permanent Filtration .....	1.1 mm Al / 75 kV IEC60522:1999
Radiation Protection (In accordance with IEC60601-1-3:2008):	
Leakage Technique Factor .....	150 kV, 3.4 mA
X-ray Coverage .....	430 × 430 mm at SID 1000 mm
Weight (Approx.) .....	24 kg
High Voltage Receptacle .....	To meet the requirements of IEC60526 Corrigendum1:2010
Cooling Method .....	Natural or forced air
Tube Housing Model Number .....	XH-112V

## Absolute Maximum and Minimum Ratings

(At any time, these values must not be exceeded.)



Maximum X-ray Tube Voltage (IEC60613:2010):	
Radiographic .....	150 kV
Fluoroscopic .....	125 kV
Between Anode (or Cathode) and Ground .....	75 kV
Minimum X-ray Tube Voltage .....	40 kV
Maximum X-ray Tube Current (IEC60613:2010):	
Large Focus .....	1000 mA
Small Focus .....	500 mA
Maximum Filament Current:	
Large Focus .....	5.8 A
Small Focus .....	5.2 A
Filament Voltage:	
Large Focus (At maximum filament current 5.8 A) .....	12.8 ~ 17.2 V
Small Focus (At maximum filament current 5.2 A) .....	10.2 ~ 13.8 V
Filament Frequency Limits .....	0 ~ 25 kHz
Continuous Anode Input Power (IEC60613:2010) .....	240 W (335 HU/s)
(Fluoroscopic, repeated radiographic or mixed exposure)	
Thermal Characteristics:	
Anode Heat Content .....	210 kJ (300 kHU)
Maximum Anode Heat Dissipation .....	710 W (1000 HU/s)
X-ray Tube Assembly Heat Content .....	1070 kJ (1508 kHU)
Nominal Continuous Input Power (IEC60613:2010):	
Without Air-circulator .....	215 W (18 kHU/min)

## Environmental Limits

Operating Limits:	
Temperature .....	10 ~ 40 °C
Humidity .....	30 ~ 85 %
(No condensation)	
Atmospheric Pressure .....	70 ~ 106 kPa
Shipping and Storage Limits:	
Temperature .....	-20 ~ 70 °C
Humidity .....	20 ~ 90 %
(No condensation)	
Atmospheric Pressure .....	50 ~ 106 kPa

## Maximum Rating Charts (Absolute Maximum Rating Charts)



Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 50Hz



Refer to IEC60613:2010

## Maximum Rating Charts (Spot-Film Rating Charts)



Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 50Hz



Refer to IEC60613:2010

## Maximum Rating Charts (Absolute Maximum Rating Charts)



Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 60Hz



Refer to IEC60613:2010

## Maximum Rating Charts (Spot-Film Rating Charts)



Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 60Hz



Refer to IEC60613:2010

## Maximum Rating Charts (Absolute Maximum Rating Charts)



Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 180Hz



Refer to IEC60613:2010

## Maximum Rating Charts (Spot-Film Rating Charts)

Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 180Hz




Refer to IEC60613:2010

## Emission & Filament Characteristics

Constant potential high-voltage generator


Nominal Focal Spot Value: 1.2 ■



Note1) For Reference Only

Note2) Refer to IEC60613:2010

Nominal Focal Spot Value: 0.6 □




Note1) For Reference Only


Note2) Refer to IEC60613:2010

## Thermal Characteristics

X-ray Tube Assembly Heating / Cooling Curve



Anode Heating / Cooling Curve




The heating curves are showing examples of average input power to the anode in operation.




### Dimensional Outline

Unit: mm




#### TERMINAL CONNECTIONS



Note) Do not connect terminal No. 1 and No. 5 or No. 6 in series circuit.

**EXPLANATION OF SYMBOLS**  
 CATHODE TERMINAL  
 C : COMMON  
 L : LARGE FOCUS  
 S : SMALL FOCUS

**TERMINAL CONNECTIONS**  
 C1 : COMMON  
 M : MAIN WINDING OF THE STATOR  
 A : AUX. WINDING OF THE STATOR  
 NC : NON-CONNECTION  
 ET : EARTH TERMINAL



↑ : CENTRAL X-RAY  
 ANODE & CATHODE TERMINAL  
 : IEC60526 TYPE

## OVERSEAS SUBSIDIARIES AND AFFILIATES



### EU REPRESENTATIVE

· **TOSHIBA ELECTRONICS EUROPE GMBH**

HANSAALLEE 181 40549 DÜSSELDORF, GERMANY  
PHONE +49 (211) 5296-107      FAX +49 (211) 5296-402

**For Sales & Technical Services, please contact the following representative:**

· **TOSHIBA ELECTRONICS EUROPE GMBH**

HANSAALLEE 181 40549 DÜSSELDORF, GERMANY  
PHONE +49 (211) 5296-107      FAX +49 (211) 5296-402

· **TOSHIBA AMERICA ELECTRONIC COMPONENTS, INC.**

2150 EAST LAKE COOK ROAD, SUITE 310  
BUFFALO GROVE, ILLINOIS 60089 USA  
PHONE +1 (847) 484-2400      FAX +1 (847) 541-7287

· **TOSHIBA ELECTRON DEVICES & MATERIALS (SHANGHAI) CO., LTD. (TEMS)**

RM1606, SH-PLAZA,  
No.336, XIZANG ROAD (MIDDLE), SHANGHAI, 200001, CHINA  
PHONE +86 (21) 6361-0077      FAX +86 (21) 6351-5760