


ROTANODE™
E7239X
E7239FX  **0197**
E7239GX

Rotating Anode X-ray Tube Assembly

- ◆ Rotating anode X-ray tube assembly for the purpose of general diagnostic X-ray procedures.
- ◆ Specially processed Rhenium-tungsten faced molybdenum target of 74 mm diameter.
- ◆ These tubes have foci 2.0 and 1.0, and are available for a maximum tube voltage 125 kV.
- ◆ Accommodated with IEC60526 type high-voltage cable receptacles.



General Data

IEC Classification (IEC60601-1:2005) Class I ME EQUIP MET

Electrical:

Circuit:

High Voltage Generator Constant Potential High-Voltage Generator
Grounding Center-grounded

Nominal X-ray Tube Voltage (IEC60613:2010):

Radiographic 125 kV

Nominal Focal Spot Value (IEC60336:2005):

Large Focus 2.0
Small Focus 1.0

Nominal Anode Input Power (at 0.1s):

	50 Hz	60 Hz
Large Focus	42.5 kW	47 kW
Small Focus	21 kW	22.5 kW

Nominal Radiographic Anode Input Power (IEC60613:2010):

	50 Hz	60 Hz
Large Focus	42.5 kW	47 kW
Small Focus	21 kW	22.5 kW

★The information contained herein is presented only as a guide for the applications of our products.
No Responsibility is assumed by TOSHIBA ELECTRON TUBES & DEVICES CO.,LTD.(TETD) for any infringements of patents or other rights of the third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of TETD or others.
★The information contained herein may be changed without prior notice. It is therefore advisable to contact TETD before proceeding with the design of equipment incorporating this product.

Motor Ratings:

Stator: XS-AV

		Starting		Running
		50/60		50/60
Driven Frequency	[Hz]	50/60		50/60
Input Power	[W]	1050	270	43
Voltage ^{1) 3)}	[V]	200	100	40
Current ²⁾	[A]	6.0	3.0	1.2
Min. Speed Up ⁴⁾	[s]	0.8	1.5	-
Capacitor	[μ F]	24	24	24

Note: 1) Applied voltage between common and main terminal.

2) Common current.

3) The every applied voltage must be never exceeded 110% of the above specification.

4) The speed-up time is allowed up to 110% of the above specification.

Anode Speed:

50 Hz	Minimum 2700 min ⁻¹
60 Hz	Minimum 3200 min ⁻¹

Stator Resistance:

Common-Main Winding	27.5 Ω
Common-Auxiliary Winding	58.0 Ω

Resistance Between Housing and Low Voltage Terminals Minimum 2 M Ω

Normal Operating Range of the Housing Temperature 16 ~ 75 °C

Mode of Operation Intermittent

Mechanical:

Dimensions See dimensional outline

Overall Length 479 mm

Maximum Diameter 152.4 mm

Target:

Anode Angle 16 degrees

Diameter 74 mm

Construction Rhenium-Tungsten faced Molybdenum

Filtration:

Permanent Filtration 0.9 mm Al / 75 kV IEC60522:1999

Available Additional Filter combination (0.4 - 1.5 mm) Maximum 2.4 mm Al / 75 kV

Radiation Protection (In accordance with IEC60601-1-3:2008):

Leakage Technique Factor 125 kV, 4 mA

X-ray Coverage 354 × 354 mm at SID 750 mm

Weight (Approx.) 16 kg

High Voltage Receptacle To meet the requirements of IEC60526 Corrigendum1:2010

Cooling Method Natural or forced air

Tube Housing Model Number:

E7239X XH-121

E7239FX XH-126

E7239GX XH-150

Absolute Maximum and Minimum Ratings

(At any time, these values must not be exceeded.)



Maximum X-ray Tube Voltage (IEC60613:2010):	
Radiographic	125 kV
Between Anode (or Cathode) and Ground	62.5 kV
Minimum X-ray Tube Voltage	40 kV
Maximum X-ray Tube Current: (IEC60613:2010)	See rating charts
Large Focus	570 mA
Small Focus	340 mA
Maximum Filament Current:	
Large Focus	5.1 A
Small Focus	5.1 A
Filament Voltage:	
Large Focus (At maximum filament current 5.1 A)	7.7 ~ 10.4 V
Small Focus (At maximum filament current 5.1 A)	5.8 ~ 7.8 V
Filament Frequency Limits	0 ~ 25 kHz
Continuous Anode Input Power (IEC60613:2010)	60 W (85HU/s)
Thermal Characteristics:	
Anode Heat Content	100 kJ (140 kHU)
Maximum Anode Heat Dissipation	475 W (667 HU/s)
X-ray Tube Assembly Heat Content	900 kJ (1250 kHU)
Nominal Continuous Input Power (IEC60613:2010):	
Without Air-circulator	180 W (15 kHU/min)

Environmental Limits

Operating Limits:	
Temperature	10 ~ 40 °C
Humidity	30 ~ 85 %
	(No condensation)
Atmospheric Pressure	70 ~ 106 kPa
Shipping and Storage Limits:	
Temperature	-20 ~ 70 °C
Humidity	20 ~ 90 %
	(No condensation)
Atmospheric Pressure	50 ~ 106 kPa



Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage
Constant Potential High-Voltage Generator
Stator Power Frequency 50Hz



Refer to IEC60613:2010

Conditions: Tube Voltage
Constant Potential High-Voltage Generator
Stator Power Frequency 60Hz




Refer to IEC60613:2010

Emission & Filament Characteristics


Constant Potential High-Voltage Generator

Nominal Focal Spot Value: 2.0 ■



Note1) For Reference Only
 Note2) Refer to IEC60613:2010


Nominal Focal Spot Value: 1.0 □




Note1) For Reference Only
 Note2) Refer to IEC60613:2010

Thermal Characteristics

X-ray Tube Assembly Heating / Cooling Curve






Anode Heating / Cooling Curve




The heating curves are showing example of average input power to anode in operation.

Dimensional Outline of E7239X

Unit mm



TERMINAL CONNECTIONS



TEMPERATURE RELAY
(NORMALLY CLOSED)

Note) Do not connect terminal No.1 and No.5 or 6 in series circuit.


EXPLANATION OF SYMBOLS

CATHODE TERMINAL

- C : COMMON
- L : LARGE FOCUS
- S : SMALL FOCUS

TERMINAL CONNECTIONS

- C1 : COMMON
- M : MAIN WINDING OF THE STATOR
- A : AUX. WINDING OF THE STATOR
- NC : NON-CONNECTION
- ET : EARTH TERMINAL




$-1.5\text{mm} \leq A \leq 1.5\text{mm}$
 $-1.5\text{mm} \leq B \leq 1.5\text{mm}$


▲ : CENTRAL X-RAY
 ANODE & CATHODE TERMINAL
 : IEC60526 TYPE

Dimensional Outline of E7239FX

Unit mm




TERMINAL CONNECTIONS



TEMPERATURE RELAY
(NORMALLY CLOSED)

Note) Do not connect terminal No.1 and No.5 or 6 in series circuit.




EXPLANATION OF SYMBOLS

CATHODE TERMINAL
 C : COMMON
 L : LARGE FOCUS
 S : SMALL FOCUS

TERMINAL CONNECTIONS

C1 : COMMON
 M : MAIN WINDING OF THE STATOR
 A : AUX. WINDING OF THE STATOR
 NC : NON-CONNECTION
 ET : EARTH TERMINAL







$-1.5\text{mm} \leq A \leq 1.5\text{mm}$
 $-1.5\text{mm} \leq B \leq 1.5\text{mm}$

▲ : CENTRAL X-RAY
 ANODE & CATHODE TERMINAL
 : IEC60526 TYPE

Dimensional Outline of E7239GX

Unit mm



Note) Do not connect terminal No.1 and No.5 or 6 in series circuit.


EXPLANATION OF SYMBOLS

CATHODE TERMINAL
 C : COMMON
 L : LARGE FOCUS
 S : SMALL FOCUS

TERMINAL CONNECTIONS
 C1 : COMMON
 M : MAIN WINDING OF THE STATOR
 A : AUX. WINDING OF THE STATOR
 NC : NON-CONNECTION
 ET : EARTH TERMINAL

▲ : CENTRAL X-RAY ANODE & CATHODE TERMINAL
 : IEC60526 TYPE

OVERSEAS SUBSIDIARIES AND AFFILIATES



EU REPRESENTATIVE

·TOSHIBA ELECTRONICS EUROPE GMBH

HANSAALLEE 181 40549 DÜSSELDORF, GERMANY

PHONE +49 (211) 5296-107 FAX +49 (211) 5296-402

For Sales & Technical Services, please contact the following representative:

·TOSHIBA ELECTRONICS EUROPE GMBH

HANSAALLEE 181 40549 DÜSSELDORF, GERMANY

PHONE +49 (211) 5296-107 FAX +49 (211) 5296-402

·TOSHIBA AMERICA ELECTRONIC COMPONENTS, INC.

2150 EAST LAKE COOK ROAD, SUITE 310

BUFFALO GROVE, ILLINOIS 60089 USA

PHONE +1 (847) 484-2400 FAX +1 (847) 541-7287

·TOSHIBA ELECTRON DEVICES & MATERIALS (SHANGHAI) CO., LTD. (TEMS)

RM1606, SH-PLAZA,

No.336, XIZANG ROAD (MIDDLE), SHANGHAI, 200001, CHINA

PHONE +86 (21) 6361-0077 FAX +86 (21) 6351-5760

 **TOSHIBA ELECTRON TUBES & DEVICES CO., LTD.**
(Product Development Engineering Department)

1385 SHIMOISHIGAMI, OTAWARA-SHI, TOCHIGI-KEN, 324-8550, JAPAN

PHONE: +81-287-26-6666 FAX: +81-287-26-6060

<http://www.toshiba-tetd.co.jp/>

·Toshiba Electron Tubes & Devices Co., Ltd. meets the Environmental Management System Standard, ISO14001.

·Toshiba Electron Tubes & Devices Co., Ltd. meets internationally recognized Standards for Quality Management System ISO9001, ISO13485.