

ROTANODE™
E7869X C E 0197

Rotating Anode X-ray Tube Assembly

- ◆ High speed rotating anode X-ray tube assembly for high energy radiographic operations.
- ◆ For the purpose of general diagnostic X-ray procedures.
- ◆ This tube has foci 1.2 and 0.6, and is available for a maximum tube voltage 150 kV.
- ◆ This tube assembly has specially processed rhenium-tungsten faced molybdenum target of 100 mm diameter anode disc and are accommodated with IEC60526 type high-voltage cable receptacles.



General Data

IEC Classification (IEC60601-1:2005) Class I ME EQUIPMENT

Electrical:

Circuit:

High Voltage Generator Constant Potential High-Voltage Generator
 Grounding Center-grounded

Nominal X-ray Tube Voltage (IEC60613:2010):

Radiographic 150 kV

Nominal Focal Spot Value (IEC60336:2005):

Large Focus 1.2
 Small Focus 0.6

Nominal Anode Input Power (at 0.1s)

See rating charts
 50 Hz 60 Hz 180 Hz

Large Focus 53 kW 58 kW 100 kW
 Small Focus 21 kW 23 kW 40 kW

Nominal Radiographic Anode Input Power (IEC60613:2010):

50 Hz 60 Hz 180 Hz

Large Focus 53 kW 58 kW 97 kW
 Small Focus 21 kW 23 kW 38 kW

★The information contained herein is presented only as a guide for the applications of our products. No Responsibility is assumed by TOSHIBA ELECTRON TUBES & DEVICES CO.,LTD.(TETD) for any infringements of patents or other rights of the third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of TETD or others.
 ★The information contained herein may be changed without prior notice. It is therefore advisable to contact TETD before proceeding with the design of equipment incorporating this product.

Motor Ratings:

Stator: XS-AG

		Starting		Running	
Driven Frequency	[Hz]	180	50/60	180	50/60
Input Power	[W]	3710	800	200	90
Voltage ^{4) 6)}	[V]	420	190	80	48
Current ⁵⁾	[A]	9.6	10.2	3.0	2.2
Min. Speed Up ^{2) 8)}	[s]	1.2	0.8	-	-
Capacitor	[μ F]	6	44	6	44
Min. Braking ^{3) 8)}	[s]	2.0 (DC 100V)			

Note 1) To be obtained with TOSHIBA starter RS-200 or equivalent.

2) The speed up time from normal speed to high speed is 2/3 times of the specified speed up time from 0 to high speed, which is described on motor rating table.

3) To be applied for high speed rotation.

4) Applied voltage between common and main terminal.

5) Common current.

6) The every applied voltage must be never exceeded 110% of the above specification.

7) No more than two high speed starts per minute are permissible.

8) The speed-up time is allowed up to 110% of the above specification.

9) The generator manufacture may choose different values. The above table is one of the recommend conditions.

Anode Speed*:

50 Hz	Minimum 2700 min ⁻¹
60 Hz	Minimum 3200 min ⁻¹
180 Hz	Minimum 9700 min ⁻¹

Stator Resistance:

Common-Main Winding	9.4 Ω
Common-Auxiliary Winding.....	28.3 Ω

Resistance between Housing and Low Voltage Terminals Minimum 2 M Ω

Normal Operating Range of the Housing Temperature 16 ~ 75 °C

Mode of Operation Intermittent

*Note: The revolution acceleration must not exceed 150 revolutions per square second.

After high speed rotation, dynamic braking must slow anode rotor to less than 3000 min⁻¹ within 10 seconds, but not sooner than 2 seconds.

Mechanical:

Dimensions	See dimensional outline
Overall Length	496 mm
Maximum Diameter	195 mm
Target:	
Anode Angle	12 degrees
Diameter	100 mm
Construction	Rhenium-Tungsten faced Molybdenum
Filtration:	
Permanent Filtration	1.1 mm Al / 75 kV IEC60522:1999
Available Additional Filter combination (0.4 - 1.5 mm)	Maximum 2.6 mm Al / 75 kV
Radiation Protection (In accordance with IEC60601-1-3:2008):	
Leakage Technique Factor	150 kV, 4 mA
X-ray Coverage	430 × 430 mm at SID 1000 mm
Weight (Approx.)	24 kg
High Voltage Receptacle	To meet the requirements of IEC60526 Corrigendum1:2010
Cooling Method	Natural or forced air
Tube Model Number	E7869
Tube Housing Model Number	XH-112V

Absolute Maximum and Minimum Ratings **(At any time, these values must not be exceeded.)**

Maximum X-ray Tube Voltage (IEC60613:2010):	
Radiographic	150 kV
Between Anode (or Cathode) and Ground	75 kV
Minimum X-ray Tube Voltage	40 kV
Maximum X-ray Tube Current (IEC60613:2010)	
Large Focus	See rating charts
Small Focus	1000 mA
Small Focus	500 mA
Maximum Filament Current:	
Large Focus	5.8 A
Small Focus	5.1 A
Filament Voltage:	
Large Focus (At maximum filament current 5.8 A)	12.4 ~ 16.7 V
Small Focus (At maximum filament current 5.1 A)	9.3 ~ 12.5 V
Filament Frequency Limits	0 ~ 25 kHz
Continuous Anode Input Power (IEC60613:2010)	300 W (420 HU/s)
Thermal Characteristics:	
Anode Heat Content	420 kJ (600 kHU)
Maximum Anode Heat Dissipation	1750 W (2465 HU/s)
X-ray Tube Assembly Heat Content	1420 kJ (2000 kHU)
Nominal Continuous Input Power (IEC60613:2010)	
Without Air-circulator	278 W (23 kHU/min)

Environmental Limits

Operating Limits:

Temperature	10 ~ 40°C
Humidity	30 ~ 85 %
	(No condensation)
Atmospheric Pressure	70 ~ 106 kPa


Shipping and Storage Limits:

Temperature	-20 ~ 70°C
Humidity	20 ~ 90 %
	(No condensation)
Atmospheric Pressure	50 ~ 106 kPa

Maximum Rating Charts (Absolute Maximum Rating Charts)


Conditions: Tube Voltage
Constant Potential High-Voltage Generator
Stator Power Frequency 50Hz

Nominal Focal Spot Value: 1.2 ■



Refer to IEC60613:2010

Nominal Focal Spot Value: 0.6 □




Refer to IEC60613:2010


Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage
Constant Potential High-Voltage Generator
Stator Power Frequency 60Hz

Nominal Focal Spot Value: 1.2 ■




Nominal Focal Spot Value: 0.6 □



Maximum Rating Charts (Absolute Maximum Rating Charts)


Conditions: Tube Voltage
Constant Potential High-Voltage Generator
Stator Power Frequency 180Hz

Nominal Focal Spot Value: 1.2 ■



Refer to IEC60613:2010

Nominal Focal Spot Value: 0.6 □




Refer to IEC60613:2010

Emission & Filament Characteristics

Constant Potential High-Voltage Generator


Nominal Focal Spot Value: 1.2 ■



Note 1) For Reference Only

Note 2) Refer to IEC60613:2010

Nominal Focal Spot Value: 0.6 □




Note 1) For Reference Only


Note 2) Refer to IEC60613:2010

Thermal Characteristics

X-ray Tube Assembly Heating / Cooling Curve





Anode Heating / Cooling Curve




Dimensional Outline


Unit: mm



TERMINAL CONNECTIONS



Note) Do not connect terminal No. 1 and No. 5 or No. 6 in series circuit.



- 1.5mm ≤ A ≤ 1.5mm
- 1.5mm ≤ B ≤ 1.5mm

EXPLANATION OF SYMBOLS

CATHODE TERMINAL

- C : COMMON
- L : LARGE FOCUS
- S : SMALL FOCUS

TERMINAL CONNECTIONS

- C1 : COMMON
- M : MAIN WINDING OF THE STATOR
- A : AUX. WINDING OF THE STATOR
- NC : NON-CONNECTION
- ET : EARTH TERMINAL

- ↑ : CENTRAL X-RAY
- ANODE & CATHODE TERMINAL : IEC60526 TYPE

OVERSEAS SUBSIDIARIES AND AFFILIATES



EU REPRESENTATIVE

·TOSHIBA ELECTRONICS EUROPE GMBH

HANSAALLEE 181 40549 DÜSSELDORF, GERMANY
PHONE +49 (211) 5296-107 FAX +49 (211) 5296-402

For Sales & Technical Services, please contact the following representative:

·TOSHIBA ELECTRONICS EUROPE GMBH

HANSAALLEE 181 40549 DÜSSELDORF, GERMANY
PHONE +49 (211) 5296-107 FAX +49 (211) 5296-402

·TOSHIBA AMERICA ELECTRONIC COMPONENTS, INC.

2150 EAST LAKE COOK ROAD, SUITE 310
BUFFALO GROVE, ILLINOIS 60089 USA
PHONE +1 (847) 484-2400 FAX +1 (847) 541-7287

·TOSHIBA ELECTRON DEVICES & MATERIALS (SHANGHAI) CO., LTD. (TEMS)

RM1606, SH-PLAZA,
No.336, XIZANG ROAD (MIDDLE), SHANGHAI, 200001, CHINA
PHONE +86 (21) 6361-0077 FAX +86 (21) 6351-5760

 **TOSHIBA ELECTRON TUBES & DEVICES CO., LTD.**
(Product Development Engineering Department)

1385 SHIMOISHIGAMI, OTAWARA-SHI, TOCHIGI-KEN, 324-8550, JAPAN
PHONE: +81-287-26-6666 FAX: +81-287-26-6060
<http://www.toshiba-tetd.co.jp/>

·Toshiba Electron Tubes & Devices Co., Ltd. meets the Environmental Management System Standard, ISO14001.

·Toshiba Electron Tubes & Devices Co., Ltd. meets internationally recognized Standards for Quality Management System ISO9001, ISO13485.